

PRECISION

Precise
Pipetting

Laboratory Exclusive


With "caring for human health"
as the initial intention.
With the mission of "improving the quality of life".

永岳医疗科技（昆山）有限公司

YONG YUE MEDICAL TECHNOLOGY (KUNSHAN)CO.,LTD



YONG YUE
永 | 岳 | 医 | 疗



COMPANY PROFILE

Yong Yue Medical Technology (Kunshan) Co., Ltd. is a manufacturing enterprise specializing in research and development, production and sales of medical devices and laboratory consumables. It has passed ISO13485 quality management system certification and European union certification 0123. Products are popular in more than 80 countries and regions.

Basic Company Data	
Established Time	2001 Year
Floor Space	19,000 m ²
Covered Area	12,000 m ²
Class 100,000 clean room	5,000 m ²
Mould Workshop	2,000 m ²
Number of staff	200

Model Display



single channel pipette



8-channel
adjustable pipette



12-channel
adjustable pipette

Good News

Yongyue can provide many stypes for Filter Tips with you.



Product Cat#	Volume	Tip Length	Description	Compatible Pipettor	
YPT0010-T-F-S	10UL	39 or 47mm	Filter Steriliza	Eppendorf Gilson Thermo Brand TianLong YongYue	
YPT0020-T-F-S	20UL	56mm			
YPT0100-T-F-S	100UL	56mm			
YPT0200-T-F-S	200UL	56mm			
YPT0300-T-F-S	300UL	59mm			
YPT1000-T-F-S	1000UL	78mm			
YPT1250-T-F-S	1250UL	100mm			
YPT5000-T-F-S	5ML	150mm			Eppendorf Thermo
YPT10000-T-F-S	10ML	160mm			

Pangaea Manual Single Channel Pipette

Classic design Lightweight and precise

- Hot-selling models in the international market, in line with color aesthetic design.
- The gun is light in weight, light and convenient to use, and can reduce hand fatigue.
- Use accessory tools for quick and easy calibration and repairs



Yongyue can provide many styles for Filter Tips with you.



Precise pipetting



Ergonomics



Strong versatility



High temperature sterilization

Finger

Pipetting is more comfortable

Tip ejector

Fast and lightweight, reducing hand fatigue

Catapult stick


Easy disassembly, chemical resistance, autoclavable

Nozzle cone

Chemical resistance, suitable for various suction cups

Hydrophobic filter element

Avoid contamination and damage inside the pipette



Spin button

Easily spin and select dispense volumes

Digital window

Clear and simple

Features:

- ✓ Ergonomic design
- ✓ Various specifications, measuring range (0.5~1000ul)
- ✓ 10ul and above are equipped with replaceable hydrophobic filter element

Yongyue can provide many styles for Filter Tips with you.

Pangaea Series Single Channel Pipettes



YONGYUE
永 | 岳 | 医 | 疗

Number	Channels	Norm	Variable component	Measure volume	Inaccuracy	Imprecision	Applicable tips
YYPT-10-1	1	0.5-10ul	0.1ul	10ul	±1.00%	0.80%	10ul
				5ul	±1.50%	1.50%	
				1ul	±2.50%	1.50%	
YYPT-20-1	1	2-20ul	0.5ul	10ul	±1.00%	0.80%	
				5ul	±1.50%	1.50%	
				1ul	±2.50%	1.50%	
YYPT-50-1	1	5-50ul	0.5ul	50ul	±0.60%	0.30%	20ul & 50ul & 100ul & 200ul
				25ul	±0.90%	0.60%	
				5ul	±2.00%	2.00%	
YYPT-100-1	1	10-100ul	1.0ul	100ul	±0.80%	0.20%	
				50ul	±1.00%	0.40%	
				10ul	±3.00%	1.00%	
YYPT-200-1	1	20-200ul	1.0ul	200ul	±0.60%	0.20%	
				100ul	±0.80%	0.30%	
				20ul	±2.50%	0.80%	
YYPT-1000-1	1	100-1000ul	5ul	1000ul	±0.60%	0.20%	1000ul& 1250ul
				500ul	±0.70%	0.25%	
				100ul	±2.00%	0.70%	
YYPT10M-1	1	1-10ml	100ul	10000ul	±0.60%	0.20%	10ml
				500ul	±1.20%	0.30%	
				1000ul	±3.00%	0.60%	

8-channel pipette for nucleic acid extraction

Sampling analysis Accurate and reliable

- High consistency between channels and good sealing with tips ensures stable results between repeated pipetting batches
- Use accessory tools for quick and easy calibration and repairs



Elastic sealing structure



Anti-slip design



Fast and light



For nucleic acid extraction



Yongyue can provide many styles for Filter Tips with you.

Finger

Pipetting is more comfortable

Digital window


Clear and simple

Dispensing head

360 degree rotation for easy pipetting

Nozzle connection

Ingenious design of nozzle connecting piece, easy to observe the sealing status of the tip



Spin button

Easily spin and select dispense volumes

Nozzle ejector

It can push out 8 or 12 tips at the same time

Features:

- ✓ 8 channels 100~1000ul suitable for 96-well plate of nucleic acid extraction instrument
- ✓ A variety of specifications can be selected, the measuring range (0.5~1000ul)

Yongyue can provide many stypes for Filter Tips with you.

Technical Parameter

Number	Channels	Norm	Variable component	Measure volume	Inaccuracy	Imprecision	Applicable tips
YYPT-10-8	8	0.5-10ul	0.1ul	10ul	± 1.50%	1.50%	10ul
				5ul	± 2.50%	2.50%	
				1ul	± 4.00%	4.00%	
YYPT-50-8	8	5-50ul	0.5ul	50ul	± 1.00%	0.50%	20ul & 50ul & 100ul & 200ul
				25ul	± 1.50%	1.00%	
				5ul	± 3.00%	2.00%	
YYPT-300-8	8	50-300ul	5.0ul	300ul	± 0.70%	0.25%	
				150ul	± 1.00%	0.50%	
YYPT-1000-8	8	100-1000ul	10.0ul	50ul	± 1.50%	0.80%	
				100ul	± 0.60%	0.20%	
				500ul	± 0.70%	0.25%	
				100ul	± 2.00%	0.70%	

Yongyue can provide many stypes for Filter Tips with you.

12-channel pipette for nucleic acid extraction

Fast and accurate Wide liquid separation range

- High consistency between channels and good sealing with tips ensures stable results between repeated pipetting batches
- Use accessory tools for quick and easy calibration and repairs



Yongyue can provide many styles for Filter Tips with you.



Strong versatility



Anti-slip design



Fast and light



For nucleic acid extraction

Finger

Pipetting is more comfortable

Digital window

Clear and simple

Dispensing head

360 degree rotation for easy pipetting

Nozzle connection

Ingenious design of nozzle connecting piece, easy to observe the sealing status



Spin button

Easily spin and select dispense volumes

Nozzle ejector

It can push out 8 or 12 tips at the same time

Features:

- ✓ Various specifications, measuring range (0.5~300ul)
- ✓ 12 channels 0.5~300ul suitable for 96-well plate of nucleic acid extraction instrument

Yongyue can provide many stypes for Filter Tips with you.

Technical Parameter

Number	Channels	Norm	Variable component	Measure volume	Inaccuracy	Imprecision	Applicable tips
YYPT-10-12	12	0.5-10ul	0.1ul	10ul	±1.50%	1.50%	10ul
				5ul	±2.50%	2.50%	
				1ul	±4.00%	4.00%	
YYPT-50-12	12	5-50ul	0.5ul	50ul	±1.00%	0.50%	20ul & 50ul & 100ul & 200ul
				25ul	±1.50%	1.00%	
				5ul	±3.00%	2.00%	
YYPT-300-12	12	50-300ul	5.0ul	300ul	±0.70%	0.25%	200ul
				150ul	±1.00%	0.50%	
				50ul	±1.50%	0.80%	

Yongyue can provide many styles for Filter Tips with you.

Laboratory
exclusive

Choose Yongyue
Yongyue will support to you

You will get ***Discount or Gift***



TEL: 400-000-9961
FAX: 0512-57260163
Email: sales@yonyue.com
ADD: No.6 Jindie Road, Kunshan City,
Kunshan 215301, PR China

URL: www.yonyue.com



永岳医疗科技（昆山）有限公司

YONG YUE MEDICAL TECHNOLOGY (KUNSHAN) CO., LTD



YONGYUE
永 | 岳 | 医 | 疗

Precise
Pipetting